

## Cartridge Chemistries (Beckman DxC,LX,CX)

Cat#	Description	Abbreviation	Method	Form	Size
L-201802	Alanine Aminotransferase Assay Kit	ALT	IFCC	Dual Liquid-stable	2×200 Tests
L-200302	Albumin Assay Kit	ALB	Bromcresol Purple Colorimetric	Mono Liquid-stable	2×300 Tests
L-202102	Alkaline Phosphatase Assay Kit	ALP	AMP	Dual Liquid-stable	2×200 Tests
L-202602	Aspartate Transaminase Assay Kit	AST	IFCC	Dual Liquid-stable	2×200 Tests
L-200102	Cholesterol Assay Kit	CHOL	Enzymatic	Dual Liquid-stable	2×300 Tests
L-201902	Creatine Kinase Assay Kit	CK	NAC	Dual Liquid-stable	2×200 Tests
L-202002	Creatine Kinase Isoenzyme Assay Kit	CKMB	NAC-Immuno-inhibition	Dual Liquid-stable	2×60 Tests
L-201102	Creatinine Assay Kit	CREA	Picric Acid	Dual Liquid-stable	2×300 Tests
L-200902	Glucose Assay Kit	GLU	Hexokinase	Dual Liquid-stable	2×300 Tests
L-202402	Lactate Dehydrogenase Assay Kit	LD	Lactate to Pyruvate	Dual Liquid-stable	2×200 Tests
L-202502	Lactate Dehydrogenase-1 Isoenzyme	LD1	NAC-Immuno-inhibition	Dual Liquid-stable	2×100 Tests
L-000203	Magnesium Assay Kit	MG	Calmagite Colorimetric	Dual Liquid-stable	2×200 Tests
L-200402	Phosphorus Assay Kit	PHOS	Ultra Violet	Dual Liquid-stable	2×300 Tests
L-200702	Total Protein Assay Kit	TP	Biuret	Mono Liquid-stable	2×300 Tests
L-200602	Triglycerides Assay Kit	TG	Enzymatic/GPO Trinder, Endpoint	Dual Liquid-stable	2×300 Tests
L-202202	Urea Assay Kit	UREA	Urease	Dual Liquid-stable	2×300 Tests
L-202302	Uric Acid Assay Kit	UA	Uricase	Dual Liquid-stable	2×300 Tests
L-201702	γ-Glutamyltransferase Assay Kit	GGT	IFCC	Dual Liquid-stable	2×200 Tests
L-201202	Apolipoprotein A1 Assay Kit	ApoA1	Immunoturbidimetric	Dual Liquid-stable	R-1: 60ml×2, R-2: 20ml×2
L-201302	Apolipoprotein B Assay Kit	ApoB	Immunoturbidimetric	Dual Liquid-stable	R-1: 60ml×2, R-2: 20ml×2
L-000104	Calcium Assay Kit	CA	Arsenazo III	Mono Liquid-stable	R: 60ml×2
L-200202	High Density Lipoprotein Cholesterol A	HDLC	Direct/	Dual Liquid-stable	R-1: 60ml×2, R-2: 20ml×2
L-200502	Lipoprotein(a) Assay Kit	Lp(a)	Immunoturbidimetric	Dual Liquid-stable	R-1: 60ml×2, R-2: 20ml×2
L-200802	Low Density Lipoprotein Cholesterol A	LDLC	Direct	Dual Liquid-stable	R-1: 60ml×2, R-2: 20ml×2
L-000704	α-Hydroxybutyrate Dehydrogenase As	HBDH	DGKC	Dual Liquid-stable	R-1: 60ml×2, R-2: 15ml×2
L-201502	α-L-Fucosidase Assay Kit	(AFU)	CNPF	Mono Liquid-stable	R: 60ml×2